

BAT™ BRIGHTNESS ACUITY TESTER



Measures glare without leaving the refracting lane.



Three variable light intensity settings simulate glare conditions.



Evaluates tinted lenses, sunglasses and anti-reflectives to quantify improvement in glare disorders.



Macular photo stress testing.



Easy to use.



Operates with rechargeable battery.



Tests all patients with any acuity chart.

BAT™ BRIGHTNESS ACUITY TESTER

SPECIFICATIONS

<i>Light Output</i>	High - 400 ft. lamberts
<i>Power Input</i>	5V DC, 750 ma.
<i>Size</i>	12" (h) x 2.4" (w) x 2.5" (d) 305mm (h) x 64mm (w) x 64mm (d)
<i>Weight</i>	1.2 lbs. (.54 kg) including battery pack

ACCESORIES

- Charger, 120 volt*
- Occluder*
- Dust Cover*
- Lamp*
- Reflector*
- Battery Pack*
- Trial Lens Holder*
- Tinted Lens Holder*
- Tinted Lens Ring Holder (10 per box)*
- Pupil Gauge Kit*



PLAY IT SAFE.

Protect your warranty. Buy your Marco product only from an authorized Marco distributor. Call Marco for the name of the distributor nearest you.



Marco-Certified Quality assures you of the very finest vision products, offering advanced technology, industry-leading durability and enhanced productivity for your practice. Marco-Certified Quality is your mark of confidence.

S047000 BR-BAT

T H E L E A D E R I N V I S I O N D I A G N O S T I C S™



11825 Central Parkway | Jacksonville, FL 32224 | 800 874-5274 | 904 642-9330 | www.marco.com